

100% Ethylene Oxide Sterilisation

Most effective sterilant  
Most efficient steriliser

AN87 / AN1087 Dosimeters

**Anprolene®**  
gas sterilisation

**EOGas®**

**Dosimeters**

Quality Assurance Products

designed for use with

Anprolene and EOGas sterilisers.

### AN87 & AN1087 Dosimeters

For use with Anprolene and EOGas cabinets running sterilisation cycles at 20°C – 30°C or at 50°C, respectively. Andersen Dosimeters provide the operative with an immediate visual assurance that the requirements for sterilisation have been fulfilled. The Dosimeters show a clear-cut colour change from orange to deep blue. Once this colour change has passed the triangle indicator on the Dosimeter, conditions for sterilising the load will have been achieved. No laboratory culture is necessary.

Andersen Dosimeters change colour in proportion to the sterilant dose which is dependent on the integration of EtO gas concentration, temperature and duration of exposure. Operatives can use this information to judge quantitatively how the conditions in the cycle have compared to those that are required to achieve sterilisation. This enables users to intelligently manage cycles to ensure that the content of the load does not comprise sterility, e.g. when sterilising large quantities of gas absorbing materials.

The Dosimeter acts as a Class 5 chemical integrator in compliance with ISO11140-1:2005

Use one Dosimeter per cycle

**Quantity:** 25 / box

**Shelf life:** 4 years from date of manufacture



*These products are not a substitute for Biological Indicators (AN2203 and AN2202).*

Anprolene, EOGas and Dosimeter are registered trademarks of H.W. Andersen Products, Inc.

**ANPRO** **ANDERSEN**  
**PRODUCTS**  
the future of gas sterilisation

H.W. Andersen Products Ltd  
Unit 808, Oakwood Business Park  
Fowler Road  
Clacton-On-Sea, CO15 4AA  
United Kingdom  
Tel.: +44 (0)1255 428 328  
email: [uk@anderseneurope.com](mailto:uk@anderseneurope.com)  
[www.anderseneurope.com](http://www.anderseneurope.com)

H.W. Andersen Products, Inc.  
Health Science Park  
3202 Caroline Drive  
Haw River, NC 27258-9564  
USA  
Tel.: 800-523-1276  
email: [customerservice@anpro.com](mailto:customerservice@anpro.com)  
Ref:- Dosimeter Flyer 201602a